

Addressing the problem of IgG quantification in biologics

M McGoldrick¹; A Murphy¹; D Kane¹;
M Hopwood¹; M McAuley¹
¹Causeway Sensors Ltd, Belfast, UK



Accurate, cost-effective and timely quantification of IgG concentration is a vital tool in the development and manufacture of antibody-based biologics. However, we are often forced to compromise when using current tools. Causeway sensors have developed a dedicated SPR based sensor platform designed to allow accurate, affordable, rapid and robust IgG quantification from a variety of sample matrices with minimal sample prep.

Wide Dynamic Range


The Standard curve of the assay operates over a range of 4-5 orders of magnitude, limiting the need for multiple dilutions and extensive sample prep to obtain an accurate result.

Accurate

Generic IgG standards often fall short of providing a universal point of accuracy in binding assay based measurement of real IgG biologics. Causeway's assay validation includes real clinical therapeutic antibodies and recombinantly expressed IgG in representative media.

Rapid


Using a 60 second per sample turnaround time allows for rapid results compatible with at-line measurement in a production environment and real-time response.




The Instrument

The Causeway instrument utilises an SPR signal with an extended linear response allowing for a wide dynamic range in measurement, granting the assay ideal performance limits for this application. The unique chips use a proprietary poly-electric layering (PEL) surface technology, demonstrating high stability and resistance to regeneration conditions, minimal non-specific interaction and is fully compatible with industry standard coupling chemistries.

Protein A-coupled chip against human monoclonal antibody



Protein G-coupled chip against rabbit monoclonal antibody



Causeway
Sensors

info@causewaysensors.com

Causeway Sensors Ltd
63 University Road
Belfast
BT7 1NF