

ONCO-Chip^{3D}: facilitating physiologically relevant drug and cell therapy screening in complex 3D tumour models

K. Paterson¹, J. Beinarovica¹, A. Sim^{1,2} and M. Zagnoni^{1,3}

¹ScreenIn3D, Glasgow, United Kingdom

²AMS Biotechnology (Europe) Ltd, Abingdon, United Kingdom

³University of Strathclyde, Glasgow, United Kingdom


contact: info@screenin3d.com

Acknowledgements

Lisa Neilson*, Emily Kay*, Joanne Edwards#, Sara Zanivan*, Chiara Braconi#

* Tumour Microenvironment and Proteomics, Cancer Research UK Beatson Institute


Institute of Cancer Science, College of Medical, Veterinary and Life Sciences, University of Glasgow, Glasgow, UK



SCREEN IN 3D

Introduction

- Cell behaviour is altered in two-dimensional (2D) cultures with respect to the native 3D *in vivo* microenvironment (Breslin & O'Driscoll 2013), which affect cellular behavior and drug responses (Tung et al. 2011; Pickl & Ries 2009).
- Complex 3D *in vitro* cancer models based on the use of patient-derived cells can facilitate immunotherapy screening and combination studies in a physiologically relevant context.
- Microfluidic technologies offer precise control over the tumor microenvironment, and require only a fraction of reagents and cells, whilst enabling an increase in throughput even when small samples are available.
- ScreenIn3D's microfluidic technology and associated protocols enable the generation, culture and screening of hundreds of 3D complex cancer models from a mixture of primary cell lines (e.g. spheroids and organoids) and biopsy tissue (e.g. tumoroids), providing at least 20X more data throughput for the same amount of material used.
- Our proprietary ONCO-Chip^{3D} technology uses a self-generated perfusion and drug concentration gradient to miniaturise the screening of limited cellular material (e.g. biopsy-derived 3D models).



Our technology

- As low as 1,000 cells in a single injection for tens of 3D experiments
- Multiple cell types can be inserted to increase model complexity at desired time points
- Shear-stress free perfusion is achieved without any external flow actuation
- 3D co-culture models, organoids and biopsy-derived fragments for several weeks-long assays
- Controllable and long-lasting molecular concentration gradients are generated to miniaturize combination studies and quantify dose response analysis

Readouts

- Cell specific response
- 3D structure markers (LGR5, LAMB1, ki67)
- Real time monitoring
- Phenotypic characterization (Growth rate, Live/dead, Immunofluorescence, Supernatant profile, Sample recovery and off-chip processing)
- Morphology, such as disaggregation or compaction (Kelm et al. 2003)
- Measure of drug efficacy using multiple parameters


Average Volume (µm³): Precise spheroid size control

Shape factor (SF): $SF = \frac{P^2}{4\pi A}$

VF: $VF = \frac{Area_{+FDA}}{Area_{-FDA}}$


Effective treatment vs Non-effective treatment

Drug concentration gradients to screen biopsy-derived tumoroids



Drug concentration gradients for combination drug screening

- Dichloroacetate (DCA) and salinomycin (SAL) were reported to exert a synergistic effect in spheroids derived from HCT116 cells ((Skeberdyte et al., 2018).
- Colorectal HCT116 cells were seeded into the device and allowed to aggregate for 24 hours.
- On day 1, drug concentration gradients were initiated and maintained for 16 hours, whilst the devices were placed in the incubator.
- 3 days later viability staining was conducted using FDA and PI.



Primary tumor fragment culture and characterization

Fragments of high grade serous ovarian tumour resections were seeded into the ONCO-Chip^{3D} and maintained in culture for 10 days, when viability staining using PI and FDA was performed.

24.08 µm, 27.52 µm, 30.96 µm, 34.40 µm, 37.84 µm, 41.28 µm, 44.72 µm

PAX8 αSMA

To demonstrate the presence of tumour cells (PAX8) and cancer-associated fibroblasts (αSMA) on day 10, the fragments were fixed in situ and stained using immunofluorescence.

Cell therapy efficacy studies

NHLF: CellTrace
CAR-T: CellTrace
MDA: unlabelled
Scale bar = 50 µm

CAR-T mediated cytotoxicity and target specificity of EGFR expressing cells in complex 3D models were assessed using combination of chemotherapy and checkpoint inhibitors, similar to clinical strategies.

amsbio University of Strathclyde Glasgow

Why use **SCREEN IN 3D** services?

- HUMAN TME
- COST EFFECTIVE
- COMBINE
- DESIGN
- PERSONALISE

References

Breslin, S. & O'Driscoll, L., 2013. Three-dimensional cell culture: the missing link in drug discovery. *Drug discovery today*, 18(5-6), pp.240-9.

Kelm, J. M., Timmins, N. E., Brown, C. J., Fussenegger, M. & Nielsen, L. K., 2003. Method for generation of homogeneous multicellular tumor spheroids applicable to a wide variety of cell types. *Biotechnol. Bioeng.* 83, 173-180.

Mulholland, T. et al. 2018. Drug screening of biopsy-derived spheroids using a self-generated microfluidic concentration gradient. *Scientific Reports*, 8 (1), 14672

Pickl, M. & Ries, C.H., 2009. Comparison of 3D and 2D tumor models reveals enhanced HER2 activation in 3D associated with an increased response to trastuzumab. *Oncogene*, 28(3), pp.461-468.

Ruppen, J. et al., 2015. Towards personalized medicine: chemosensitivity assays of patient lung cancer cell spheroids in a perfused microfluidic platform. *Lab Chip*, 15(14), pp.3076-3085.

Skeberdyte, A., Sarapinen, I., Aleksander-Krasko, J., Stankevicius, V., Sužiedelis, K., & Jarmalaitė, S. (2018). Dichloroacetate and Salinomycin Exert a Synergistic Cytotoxic Effect in Colorectal Cancer Cell Lines. *Scientific Reports*, 8(1), 1-13. <https://doi.org/10.1038/s41598-018-35815-4>

Tung, Y.-C. et al., 2011. High-throughput 3D spheroid culture and drug testing using a 384 hanging drop array. *The Analyst*, 136(3), pp.473-478.