

Combinatorial Drug Screening On 3D Ewing Sarcoma Spheroids Using Droplet-Based Microfluidics

Raphaël F.-X. TOMASI¹, Gaëtan MARY¹, Romain FEVRE²⁻³, Nadia VERTTI-QUINTERO³, Aude DURAND², Elaine DEL NERY⁴, Charles N. BAROUD²⁻³

¹: Okomera, iPEPS, the HealthTech Hub, Paris Brain Institute, Hôpital Pitié Salpêtrière, 47, Boulevard de l'Hôpital, 75013 Paris, France

²: LadHyX, CNRS, Ecole Polytechnique, Institut Polytechnique de Paris, 91120 Palaiseau, France

³: Institut Pasteur, Université Paris Cité, Physical microfluidics and Bioengineering, 25-28 Rue du Dr Roux, 75015 Paris, France

⁴: Biophenics High-Content Screening Laboratory, Translational Research Department, PICT-IBISA, Institut Curie, PSL Research University, Paris, France




Abstract Okomera develops a unique droplet microfluidic workflow for miniaturized, automated and functional cell assays in precision cancer medicine. Here we applied concomitant and sequential drug combination of two chemotherapies, etoposide and cisplatin, to the Ewing sarcoma cell line A673. Synergistic drug combination was found, notably when etoposide was applied 24 h before cisplatin.

Droplet 3D cell culture and testing





Single Drug Toxicity



Similar IC50 measured in ultra low attachment multiwell plates

Combinatorial drug screening



Funding



Contact: communications@okomera.com