

# BRAIN ORGANOID-ON-CHIP PLATFORM FOR NEUROTOXICITY EVALUATIONS




H. Castiglione<sup>1,2,3</sup>, J. Rontard<sup>3</sup>, L. Madrange<sup>1,2</sup>, T. Lemonnier<sup>1,2</sup>, C. Baquerre<sup>3</sup>, B.G.C. Maisonneuve<sup>3</sup>, T. Honegger<sup>3</sup>, F. Yates<sup>1,2</sup>, P.-A. Vigneron<sup>1,2</sup>

<sup>1</sup>SupBiotech, Ecole d'ingénieurs en Biotechnologies, Villejuif, France; <sup>2</sup>CEA/IBFJ/SEPIA, Fontenay-aux-Roses, France; <sup>3</sup>NETRI, Lyon, France; heloise.castiglione@netri.com; pierre-antoine.vigneron@supbiotech.fr

## Cerebral organoid



## Microfluidic device (NETRI)




## **Brain-Organoid-on-Chip**

- Enhanced cerebral organoid reproducibility
- Promising technology to facilitate organoid scalability & to overcome current limitations due to high organoid heterogeneity [1]

## **Neurotoxicity evaluations**

- Establishment of Scorings for quality control & toxicological evaluation of compound-exposed cortical organoids
- Development of a **Prediction Algorithm** for compound neurotoxicity evaluation

## MATERIALS AND METHODS



#### Adapted to 3D cell culture:

Two compartments separated by a porous membrane:

- Open well for 3D culture
- Perfusion channel

Adapted to industrial transfer

(5x6x4.6 mm) Porous > membrane Perfusion

**Vanillin** (CAS n°121-33-5): 100, 1 000, 10 000 nM **Biphenyl-2-ylamine** (CAS n°90-41-5): 20, 200, 2 000 μM (meOH)

**Controls:** non-exposed and vehicle-exposed organoids

200 μM (meOH)

## **QUALITY SCORING**

- Cortical organoid characterization at D+60
- Scoring scale: 5 to 0 (from most to least optimal)




### **EXPOSURE SCORING**

- For compound-exposed cortical organoids (acute & chronic)
- Compared to controls



#### PREDICTION ALGORITHM


For compound classification into 3 neurotoxicological categories



**Potential Concern** 

Compound


• Presence of a **cyst** and a large zone without cells



progenitors (SOX2), neurons (TUBB3, MAP2), and astrocytes (GFAP) (Thunder microscope,

### Higher apoptosis & DNA damage levels compared to controls:





Immunofluorescence staining of apoptosis (CCASP3) and DNA damage (yH2AX) (Thunder microscope, Leica, objective 20X).

## microscope, Leica, objective 20X).

CONCLUSION

and astrocytes (GFAP) (Thunder microscope,

Similar apoptosis & DNA damage

Immunofluorescence staining of apoptosis

(CCASP3) and DNA damage (yH2AX) (Thunder

D+61

Leica, objective 20X, circle: rosette).

Control

levels as controls:

- Brain Organoid-on-Chip platform + Scorings + Prediction Algorithm: adapted to neurotoxicity evaluations
- Vanillin exposures: no discernable impact on morphology, cytoarchitectures & viability >> low concern

Cytotoxicity

• Biphenyl-2-ylamine exposures: altered morphology & disorganized cytoarchitectures in a dose-response manner → high concern •

at least similar to controls

markers similar to controls

**Low Concern** 

Compound

Yes

Apoptosis and DNA damage ! At least one score < 4

All scores ≥ 4

**PERSPECTIVES** 

**High Concern** 

Compound

- Implementation of additional criteria organoid for cytotoxicity characterization
- Paves the way for neurotoxicological studies & drug screening

